



Puros® Block Allograft



Changing the concept
of block grafting.

1 DOCUMENTED SUCCESSFUL BONE REGENERATION

- Documented graft and implant success rates¹⁻³ make it a viable alternative to autogenous block grafting⁴⁻⁵
- Multicenter data reported 91% *Puros* Block Allograft and 99% *Tapered Screw-Vent*® Implant survival at three years¹
- Cortico-Cancellous structure retains the remodeling capabilities of cancellous bone with the strength advantages of cortical bone needed for dimensional ridge augmentation¹⁻³

2 TIME-SAVING

- Saves time, helps to reduce pain, and can shorten the patient's rehabilitation period by eliminating the need to harvest an autogenous graft

3 TUTOPLAST® PROCESS

- Sterilized and preserved using the proprietary *Tutoplast* Process, *Puros* Block Allograft offers a high-quality option for successful bone regeneration
- Retains osteoconductive properties due to the preservation of the natural bone matrix collagen and mineral composition, trabecular pattern, original porosity^{1-2,6}

Puros Allografts – Filling nature's void.
A comprehensive line of allografts for bone and soft tissue augmentation needs.



zimmer | dental


A clinically proven documented solution for restoring volume to severely resorbed ridges.¹⁻³

The unique Tutoplast Process


The proprietary *Tutoplast* Process assures the highest standard of tissue safety and quality with minimal risk of disease transmission.^{6,8}

The process preserves the valuable collagen matrix and tissue integrity while inactivating pathogens and gently removing unwanted materials, such as cells, antigens and viruses.^{6,8} The result is safe, biocompatible tissue.

For over 40 years, *Tutoplast* processed tissues have been safely used in more than three million procedures.⁸




Delipidization




Osmotic treatment



Oxidative treatment



Solvent dehydration



Low-dose gamma irradiation

Take a closer look



Figure A Preoperative: narrow ridge.



Figure B Puros Block Allograft in place.



Figure C Postoperative: Puros Block Allograft incorporated in bone after healing.



Figure D Rehabilitated ridge with implant restoration.

Clinical photos © 2012 J. Daulton Keith, Jr., DDS, FICD. All rights reserved. Individual results may vary.

Ordering information

Catalog Number	Description
68220	Puros Block Allograft, 10mm
68221	Puros Block Allograft, 15mm

To learn more about **Puros Block Allograft**, please visit us online at www.zimmerdental.com or to speak to a sales representative, call 1 (800) 854-7019.

¹ Keith JD Jr, Petrungaro P, Leonetti JA, Elwell CW, Zeren KJ, Caputo C et al. Clinical and histologic evaluation of a mineralized block allograft: results from the developmental period (2001-2004). *Int J Periodontics Restorative Dent*. 2006;26:321-327.

² Leonetti JA, Koup R. Localized maxillary ridge augmentation with a block allograft for dental implant placement: case reports. *Implant Dent*. 2003;12:217-226.

³ Keith JD Jr, Salama MA. Ridge preservation and augmentation using regenerative materials to enhance implant predictability and esthetics. *Compend Contin Educ Dent*. 2007;28:614-623.

⁴ Schwartz-Arad D, Levin L, Sigal L. Surgical success of intraoral autogenous block onlay bone grafting for alveolar ridge augmentation. *Implant Dent*. 2005 14:131-138.

⁵ Levin L, Nitzan D, Schwartz-Arad D. Success of dental implants placed in intraoral block bone grafts. *J Periodontol*. 2007;78:18-21.

⁶ Schoepf C. Allograft safety: efficacy of the Tutoplast® process. *International Magazine of Oral Implantology*. 2006;1:10-15.

⁷ von Arx T, Buser D. Horizontal ridge augmentation using autogenous block grafts and the guided bone regeneration technique with collagen membranes: a clinical study with 42 patients. *Clin Oral Impl Res*. 2006;17:359-366.

⁸ Data on file with RTI Biologics, Inc.

To receive our eNews visit us at http://www.zimmerdental.com/news_eNewsletterSignUp.aspx

For more information about our Products, Regenerative Materials and Educational Opportunities, contact us:

In the U.S. 1 (800) 854-7019
To fax an order 1 (888) 225-2483
Outside the U.S. +1 (760) 929-4300
Australia +61 (0)2 9950 5434 or 1 (800) 241 916
Canada +1 (905) 567-2073 or 1 (800) 265-0968
China +86 21 22115147

France +33 (0)1 45 12 35 35
Germany +49 (0)761 1 56 47 0
Israel +972 (0)3 6124242
Italy +39 043 855 5573
Spain +34 93 846 05 43



Zimmer Dental
1900 Aston Avenue
Carlsbad, CA 92008-7308
USA

www.zimmerdental.com