



CopiOs Extend™ Membrane

Synergy of Space and Time

1 EXTENDED RESORPTION PROFILE


- Lasts 6+ months,¹ designed to allow implant placement while providing ample time for regeneration

2 SPACE MAINTAINING CONFORMANCE

- Designed to conform to the defect with enough structural integrity for space creation

3 VERSATILE

- Demonstrated efficacy in GBR procedures with and without primary closure²



Long lasting, resorbable collagen membrane for GBR and GTR procedures

Copios Extend™ Membrane

- **Biocompatible:** Highly purified, intact porcine dermis¹
- **Easy to Use:** Supplied sterile and implantable dry or briefly hydrated
- **Convenient Handling:** Conformable and easy to reposition in the defect



INDICATIONS

Copios Extend Membrane is intended for use in oral surgical procedures as a resorbable material for:

- Use in augmentation around implants placed in immediate and delayed extraction socket
- Localized ridge augmentation for later implantation
- Alveolar ridge reconstruction for prosthetic treatment
- Filling of bone defects
- Guided bone regeneration in dehiscence defects
- Guided tissue regeneration procedures in periodontal defects

Take a closer look



Figure A Clinical image showing the extensive bone loss on the palatal area extending past apex.



Figure B *Copios Extend* was laid over the graft and tacked under the flap.



Figure C Flap sutured over the membrane.



Figure D Two months postoperative. Excellent healing.

Clinical photos ©2013 Dr. El Chaar
All rights reserved. Individual results may vary.

To learn more about the *Copios Extend* Membrane, please visit us online at www.zimmerdental.com, or to speak to a sales representative call 1 (800) 854-7019 (in the US) or +1 (760) 929-4300 (outside the US).

Ordering information

Catalog No.	Description
0190Z	<i>Copios Extend</i> Membrane, 15 x 20 mm
0191Z	<i>Copios Extend</i> Membrane, 20 x 30 mm
0192Z	<i>Copios Extend</i> Membrane, 30 x 40 mm

1. Data on file with Collagen Matrix, Inc.
2. El-Chaar ES, Wen HB. Clinical evaluation of a slow-resorbing porcine dermis membrane for implant site development: a pilot case series with and without primary closure. Presented at the Annual Meeting of the American Academy of Periodontology, September 19-22, 2014, San Francisco, CA.
©2014 Zimmer Dental Inc. All rights reserved. ZD1134, Rev. 9/14. *Copios Extend* Membrane is manufactured by Collagen Matrix, Inc. and distributed by Zimmer Dental Inc. Please note that not all products and regenerative materials are registered or available in every country/region. Please check with a Zimmer Dental representative for availability and additional information.
This document is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part. Please refer to the package inserts for important product information, including, but not limited to, contraindications, warnings, precautions, and adverse effects.

To receive our eNews visit us at http://www.zimmerdental.com/news_eNewsletterSignUp.aspx

For more information about our Products, Regenerative Materials and Educational Opportunities, contact us:

In the U.S. 1 (800) 854-7019
To fax an order 1 (888) 225-2483
Outside the U.S. +1 (760) 929-4300
Australia +61 (0)2 9950 5434 or 1 (800) 241 916
Canada + 1 (905) 567-2073 or 1 (800) 265-0968
Chile +562 231 5185

China +86 21 22115147
France +33 (0)1 45 12 35 35
Germany +49 (0)761 1 56 47 0
Israel +972 (0)3 6124242
Italy +39 0438 37681
Spain +34 93 846 05 43



Zimmer Dental
1900 Aston Avenue
Carlsbad, CA 92008-7308
USA

www.zimmerdental.com